



GelView4000 Lite

Smart Vision Gel Imaging System

Biolight Biotechnology Co., Ltd

In partnership with Verulam Scientific Ltd

Born For Gel Imaging

Features



High resolution camera, clearer imaging, capturing every detail



Breathing logo, indicating operational status of the device



A pebble shape, making it the most beautiful item in the lab



Only the size of two standard basketballs, greatly saves the laboratory space



Gesture control: Wave hand to open or close the lid, used to control the status of the lid

Voice control: Support using "Hi,Blueten" to wake up the device and other voice control operations

Device control: After connect PC, Software can be used to control imaging



Device Safety: equipped with temperature control module, the light source will close automatically when detecting overheat, avoiding high temperature damage to the device

Sample Safety: the system can set the UV light source to shut down at a certain time, avoiding sample degradation after long-term irradiation



▲ Intelligent Operation Analysis System

Product Parameters

Camera	High resolution low light digital camera	Light Source	36 LED lamp beads with a wavelength of 302nm
Resolution	3072×2048	Filters	590nm
Dynamic Range	>3.4log	Size	310 mm× 520 mm×270 mm (W×D×H)
Imaging View	13cm×13cm	Weight	4.8kg

